



A typical setup for a PET scanner application using a PhotoniQ, a SIB2316 and two silicon photomultiplier arrays (SPMArray) is shown below. The SPMArrays are positioned to detect incoming light from a scintillator crystal or optical assembly and connected to the SIB2316 by two FPC cables. The sensor interface board cable (SIB cable) connects the 32 detector outputs from the SIB2316 to a PhotoniQ IQSP480 or IQSP580 32 channel data acquisition system. Bias to the detectors is provided from bias supplies on the SIB2316 and controlled through the PhotoniQ graphical user interface. The two independent bias supplies on the SIB2316 allow each SPMArray to be individually trimmed to optimize gain matching of the silicon photomultiplier devices. The trigger output from the SIB2316 supplies the trigger to the PhotoniQ when coincident pulses exceeding a user-programmed energy threshold are detected on the two SPMArrays. Alternatively, the coincidence function can be bypassed altogether and the PhotoniQ triggered when a single pulse from either SPMArray exceeds the energy threshold. The energy threshold is also controlled through the PhotoniQ GUI. Each time the PhotoniQ is triggered, the total energy from each channel from each SPMArray (32 channels total) is digitized in parallel and sent to a PC over a USB 2.0 connection for display, logging, or real time processing.



In a PET application, the PhotoniQ is configured in pre-trigger mode so that the unit can capture the energy information from the detectors that occurred just prior to the trigger signal. This powerful triggering mode allows the data acquisition unit to be timed to the pulse peaks yet still collect the entire signal from the particle event — including the portion of the signal that preceded the peak of the event. Timing for this mode is shown below.





Vertilon Corporation has made every attempt to ensure that the information in this document is accurate and complete. Vertilon assumes no liability for errors or for any incidental, consequential, indirect, or special damages including, without limitation, loss of use, loss or alteration of data, delays, lost profits or savings, arising from the use of this document or the product which it accompanies.

Vertilon reserves the right to change this product without prior notice. No responsibility is assumed by Vertilon for any infringements of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under the patent and proprietary information rights of Vertilon Corporation.

© 2011 Vertilon Corporation, ALL RIGHTS RESERVED

AN3321.1.2 Apr 2011

Vertilon Corporation, 66 Tadmuck Road, Westford, MA 01886 / Tel: (978) 692-7070 / Fax: (978) 692-7010 / www.vertilon.com